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Tab.1 The retention time, qualitative and quantitative ion, linear regression equation and correlation coefficient of VOCs

compounds	retention time	qualitative ion m/z	quantitative ion m/z	linear regression equation	RSD%
Benzene (苯)	8.991	78	78	$Y=5.03e^{-001}X$	4.2%
Toluene (甲苯)	11.492	91	91	$Y=6.73e^{-001}X$	2.5%
p/m- Xylene (对、间二甲苯)	13.454	91	91	$Y=1.39e^{+000}X$	3.0%
o - Xylene (邻二甲苯)	13.884	91	91	$Y=8.12e^{-001}X$	3.1%
Ethylbenzene (乙苯)	13.326	91	91	$Y=8.34e^{-001}X$	7.3%
Styrene (苯乙烯)	13.834	91	91	$Y=4.77e^{-001}X$	4.7%

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Tab.2 The classification of automobile interior decorative materials

classification	decorative materials
rubber (橡胶)	sealer (车门、车窗密封胶条)
leather (皮革)	coated parts (座椅表皮、皮内饰)
fabrics (纤维纺织品)	carpet(地毯)、automobile roof(汽车顶篷)、textile (化纤纺织品或合成纤维无纺布)
silencer (消音材料)	Sponge (海绵、泡绵)
plastic (热塑性材料)	instrument panel (仪表盘)
adhesive (粘合剂)	
wooden hardboard (板材)	wood chip and slab (轿车车厢的隔板、门内装饰板)

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[1] VDA 278. Thermodesorption analysis of Volatile Organic Compounds emissions for the characterization of non-metallic automobile interior materials [S]. German Association of the Automotive Industry Recommendation , 2002

[2] VDA 277. Determination of emissions of Volatile Organic Compounds of non-metallic materials from vehicles interior [S]. German Association of the Automotive Industry Recommendation , 1995.

[3] Guangzhou Honda Automobile Co. Ltd., Specification test method for VOC [S], Guangzhou Honda Automobile Co. Ltd. Standard, DW G NO. 0094Z-SNA-000.

[4] SUZUKI Motor Corporation, Standard for control of Volatile Organic Compounds emission [S], SUZUKI Engineer Standard, SES N 2403, 2005

[5] MAZDA Motor Corporation, Standard for Volatile Organic Compounds (VOC) from parts in vehicle compartment [S], MAZDA Engineer Standard, MES CF 080, 2004

[6] JIS A 1901, Small chamber method – determination of emission of Volatile Organic Compounds and Aldehyde for building products [S], Japanese Industrial Standards, 2003

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China State Bureau of Standards, GB 18580-2001 Indoor decorating and refurbishing materials-Limit of formaldehyde emission of woodbased panels and finishing products [S], Beijing: Standards Press of China, 2006 (in Chinese)

[8] EPA TO15. Compendium of methods for the determination of toxic organic compounds in ambient air. Second Edition. Compendium Method TO-15 Determination of Volatile Organic Compounds (VOCs) in air collected in specially-prepared canisters and analyzed by Gas Chromatography-Mass Spectrometry (GC/MS), 1999

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Wu Jie-shan, Liu Chong-hua, Zhang Xiao-li, et al. Determination of volatile alkyl benzenes in automobile interior decoration materials by thermodesorption Gas Chromatography-Mass Spectrometry[J]. Journal of Instrumental Analysis. 2008, 27 (8): 874-877 (in Chinese)

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Zeng Ming, Lin Zhi-fen, Yin Da-qiang, et al. Progress on joint effect of mixture pollutants[J]. Environmental Science & Technology, 2009, 32 (2) : 86-92 (in Chinese)

**Determination and Analysis of Volatile Organic Compounds in Automobile Interior
Decoration Materials by GC-MS**

He Xiao-feng¹ Liu Yan-lin¹ Wang Gui-xia² Zhong Wen-jun²

(1 Shenzhen Institute of Information Technology, Shenzhen 518029, P. R. China)

(2 Shenzhen GD Air Testing, Inc , Shenzhen 518057, P. R. China)

Abstract: The volatile organic compounds releasing from the automobile interior decoration materials were collected by a constant temperature chamber. The VOC's concentration of the samples was pre-enriched for three times by Nutech Pre-enrichment Instrument. Then the volatile organic compounds such as benzene, toluene, xylene, ethylbenzene, styrene, tetradecane were measured by GC-MS (gas chromatography - mass spectrometry). The detection limits of this method are less than 2.0 $\mu\text{g}/\text{m}^3$. Five sorts of automobile interior decoration materials such as leather coated part, rubber sealer, textile, sponge and instrument panel were chosen for the determination of volatile organic compounds, and the results were analyzed and discussed. The results showed that five sorts of automobile interior decoration materials released the volatile organic compounds to varying degrees. The instrument panel released the most TVOC, the rubber sealer released more TVOC, the leather coated part and textile released more toluene and the sponge released more xylene. To explore and establish effective detection methods of automobile interior decoration materials is very meaningful for the choice of low toxicity and environmentally friendly automobile interior decoration materials.

Keywords: automobile interior decoration materials; volatile organic compounds; GC-MS (gas chromatography - mass spectrometry)

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